

# Neonatal Incubator Parts Drying



**Drying incubator components is an integral part of the infection control process when performing the routine cleaning and disinfection of neonate incubators. If parts aren't thoroughly dry when the incubator is put back together there is a significant risk of re-growth of harmful bacteria within the incubator - the last thing anybody wants where pre-term babies are involved.**

- Temperature range: 40°C to 90°C
- Repeatable run conditions
- 3 year warranty\*
- Standard and bespoke sizes



The range includes bench top and floor standing models as well as bespoke designed ovens like the one shown here on wheels. Moist air is expelled from the rear of the dryer.

These ovens are manufactured with a powder coated exterior and a medical grade stainless steel interior. A stainless steel drip tray is also fitted. Three heavy duty shelves are supplied with each dryer.

*This model can dry a complete neonatal incubator at one go. Ref: HDH-NEO-DRY-760-SS-ECO-WH*

The control system uses a PID digital display controller that shows both set and measured drying temperatures to give staff reassurance that the drying process is working correctly.

An over temperature safety cut out is fitted to protect expensive incubator parts from heat damage and our unique air flow system keeps the drying temperature stable and uniform at all times.

## Standard Models

<u>2 doors</u>	<u>Shelves</u>	<u>Int. HxWxD mm</u>	<u>Ext. HxWxD mm</u>
HDH-NEO-DRY-425-SS-ECO	3	580 x 960 x 690	1005 x 1065 x 800
HDH-NEO-DRY-500-SS-ECO	3	695 x 960 x 690	1120 x 1065 x 800
HDH-NEO-DRY-750-SS-ECO-WH	3	970 x 960 x 770	1495 x 1065 x 880
HDH-NEO-DRY-760-SS-ECO-WH	3	970 x 1020 x 770	1495 x 1075 x 850

<u>1 door</u>	<u>Shelves</u>	<u>Int. HxWxD mm</u>	<u>Ext. HxWxD mm</u>
HDV-NEO-DRY-350-SS-ECO	3	740 x 690 x 650	1165 x 795 x 760
HDV-NEO-DRY-500-SS-ECO-WH	3	1095 x 690 x 650	1520 x 795 x 760

*Options : moisture extraction unit, audible overheat warning, timers & stands with or without wheels*

*Bespoke sizes are available - please call us. Electrics: 230v, 50Hz, 3Kw (standard UK plug)*

**Other products include: resuscitaires, neonate incubators, radiant warmers and neonate phototherapy systems**

**Contact us on: 01663-735494**

**Email: [info@qedscientific.com](mailto:info@qedscientific.com)**