



BREAST MILK REFRIGERATORS

OPTIMUM STORAGE OF BREAST MILK

DIFFERENT SIZES, SAME QUALITY

Our breast milk refrigerators are specifically designed to store human milk in optimum conditions.

The fridge defrosts automatically and is fitted with a 2m plug for ease of installation. A temperature range of 1-5°C, but set at 3°C, and stabilised by a high-quality temperature sensor, which keeps the internal temperature even and stable.

The fridge records minimum and maximum temperatures when plugged into the power source but will continue to operate if this fails by running on the batteries installed in the panel on the side of the cabinet.

The interior of the cabinet is painted aluminium and there are plastic-coated wire shelves for ease of cleaning.

The LED display controller can be easily set to your required temperature.

A bright LED strip lamp gives clear illumination to the interior.

If the temperature drops below or raises above the set level, an audible alarm will sound.

Similarly, the alarm will sound if the door is accidentally left open or not shut properly.

The refrigerator has adjustable feet to ensure stability and castors are fitted to the rear of the product for ease of installation or relocation if required.

FEATURES

- iCool built-in temperature probe accurately monitors and stabilises drug temperatures
- External digital controller and temperature display with battery back-up
- 24/7 automatic temperature data logging with min/max temperature recording and display
- Audible and visual temperature alarms; power failure and door open alarms
- LED energy-efficient lighting
- Door lock (2 keys)
- R600a environmentally-friendly refrigerant
- 2 years warranty

Contact us to discuss your requirements and specification or to request further information.

Office

info@qedscientific.com

01663 735494

OUR RANGE

Available in three different sizes, we have a solution for your requirements.
Contact us to discuss your needs

| MODEL | LITRES | SHELVES | EXTERNAL (MM) H X W X D | WEIGHT EMPTY (KG) |
|---------|--------|---------|----------------------------|-------------------|
| BMR-60 | 60 | 3 | 725 x 445 x 485 | 28 |
| BMR-150 | 150 | 3 | 820 x 595 x 595 | 38 |
| BMR-300 | 300 | 6 | 1500 x 595 x 650 | 64 |

BMR60

- Auto defrost
- Internal fan
- Range 1°C - 5°C
- 2 + 1 small shelf
- 725h x 445w x 485d (mm)



BMR150

- Auto defrost
- Internal fan
- Range 1°C - 5°C
- 2 + 1 small shelf
- 820h x 595w x 595d (mm)

BMR300

- Auto defrost
- Internal fan
- Range 1°C - 5°C
- 6 shelves
- 1500h x 595w x 650d (mm)



