



MEDICAL WARMING CABINETS

FLUID, BLANKET AND CONTRAST MEDIA

UK MANUFACTURED CABINETS

QED Scientific are proud to offer a wide range of cabinets for warming fluids in operating theatres and critical care areas, blankets in A&E and theatre recovery and contrast media in X-Ray CT scanner rooms.

All have easy to use advanced controls such as digital temperature control with automatic setting of overheat safety and an independent protection feature.

The chambers are medical grade polished stainless steel for easy cleaning to aid infection control procedures.

Factory settings for maximum temperature are:

- Fluid warming cabinets** 60°C
- Blanket warming cabinets** 80°C
- Media warming cabinets** 39°C

The internal fans ensure fast temperature recovery after door openings.

- AWARD WINNING** - Recognised market leading design with adaptive setting of protection level
- WIDE RANGE** - 15 fluid models, 10 blanket models, 13 combined models, 7 media warming models
- EASY TO USE** - Delivered as 'Plug and play', Automatic safety protection
- RELIABLE** - 5 year warranty as standard
- 5 YEAR UK WARRANTY** - Year 1 - parts and labour, year 2-4 - parts



FLUID WARMERS

Designed to hold bottles from 1 litre up to 3 litres, or bags from half a litre to 3 litres. Capacities from 9 to 240 1 litre bottles on heavy duty wire shelves; 2 to 60 3 litre bags in wire baskets. These have a temperature range from above ambient to 60°C. Accessories include stands for for under bench models, lockable wheels for larger cabinets and door locks to prevent unauthorised access.

BLANKET WARMERS

For warming cotton blankets in theatre recovery, A&E resus rooms and critical care areas. Cabinets are supplied with baskets for storing from 4 to 50 standard cotton blankets. The temperature for warming cotton blankets is limited to 80°C. Accessories include stands or wheels for smaller models, door locks to prevent unauthorised access

COMBINATION WARMERS

For smaller operating theatres or A&E units where warm fluids and blankets are required but floor space is limited. These are based on a combination of a fluid warmer stacked on top of a blanket warmer And fixed together. Lockable wheels are fitted to all models. Common accessories include digital door locks.

CONTRAST MEDIA WARMERS

For warming bottles of contrast media in CT radiography rooms. Models are available to warm from 10 to 120 50ml bottles stored on perforated shelves. Bottles up to 200ml in size can be warmed and stored. Door locks are often requested to be fitted to these cabinets. Accessories also include wall fixing brackets and wall shelves.

Contact us to discuss your requirements and specification or to request further information and range on any of our warming cabinets.

info@qedscientific.com
01663 735494

